

[KO Validated] Anti-Sirt1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q96Eb6

Item Number: CM0021878



Introduction

Rabbit polyclonal antibody against SIRT1 (NAD-dependent deacetylase sirtuin-1), validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Suitable for studying metabolism, apoptosis, chromatin regulation, and aging.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SIRT1 Rabbit pAb
Remarks/Alias	SIR2; SIR2L1; SIR2alpha; SIRT1
Species	Human
GeneID (Human)	23411
GeneID	23411
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96EB6
Protein Weight	82kDa
Shipping	Ice bag