

[KO Validated] Anti-Sirt3 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Ntg7

Item Number: CM0002874



Introduction

SIRT3 rabbit monoclonal antibody validated for Western blot and ELISA applications. Targets human mitochondrial deacetylase SIRT3, involved in metabolism and stress response. Uniprot Q9NTG7.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SIRT3 Rabbit mAb
Remarks/Alias	SIR2L3; T3
Species	Human
GeneID (Human)	23410
GeneID	23410
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 300-399 of human SIRT3 (NP_036371.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NTG7
Protein Weight	44kDa
Shipping	Ice bag