

Ko Validated Anti-Slc25A22 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H936

Item Number: CM0019832

Introduction

KO validated rabbit polyclonal antibody targeting mitochondrial glutamate carrier SLC25A22. Validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Suitable for neuroscience, diabetes, and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SLC25A22 Rabbit pAb
Remarks/Alias	GC1; DEE3; GC-1; EIEE3; NET44; 22
Species	Human
GeneID (human)	79751
GeneID	79751
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 128-188 of human SLC25A22 (NP_078974.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H936
Protein Weight	34kDa
Shipping	Ice bag