

[Ko Validated] Anti-Slc25A23 Polyclonal Antibody For Wb, Elisa - Q9Bv35

Item Number: CM0019467

Introduction

Rabbit polyclonal antibody against SLC25A23, validated for WB and ELISA. SLC25A23 is a mitochondrial adenyl nucleotide antiporter, suitable for energy metabolism and calcium signaling research. Reacts with human, mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SLC25A23 Rabbit pAb
Remarks/Alias	APC2; MCSC2; SCAMC3; SCaMC-3; 23
Species	Human
Gene ID	79085
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human SLC25A23 (NP_077008.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BV35
Protein Weight	52kDa
Shipping	Ice bag