

[Ko Validated] Anti-Slc25A4/Ant1 Rabbit Polyclonal Antibody - P12235

Item Number: CM0020749

Introduction

KO-validated rabbit polyclonal antibody targeting SLC25A4/ANT1, the mitochondrial ADP/ATP translocase. Validated for WB, IF/ICC, ELISA; reacts with human, mouse, rat. Molecular weight ~33 kDa.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SLC25A4/ANT1 Rabbit pAb
Remarks/Alias	T1; ANT; AAC1; ANT1; PEO2; PEO3; ANT 1; PEOA2; MTDP512; MTDP512A
Species	Human
GeneID (Human)	291
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 101-200 of human SLC25A4/ANT1 (NP_001142.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12235
Protein Weight	33kDa
Shipping	Ice bag