


Ko Validated Anti-Smad1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q15797

Item Number: CM0023723



Introduction

KO-validated anti-Smad1 rabbit polyclonal antibody validated for Western blot, IHC-P, IF/ICC, and ELISA. Recognizes human, mouse, and rat Smad1. Suitable for BMP signaling, embryonic development, and cell differentiation research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Smad1 Rabbit pAb
Remarks/Alias	BSP1; JV41; BSP-1; JV4-1; MADH1; MADR1; d1
Species	Human
Gene ID (Human)	4086
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-240 of human Smad1 (NP_001003688.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15797
Protein Weight	52kDa
Shipping	Ice bag