

[Ko Validated] Anti-Smad3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Chip, Elisa - P84022

Item Number: CM0002366

Introduction

KO validated rabbit monoclonal antibody against Smad3 (P84022). Detects human, mouse, and rat Smad3. Suitable for WB, IHC-P, IP, ChIP, and ELISA. Ideal for TGF-beta signaling and transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Smad3 Rabbit mAb
Remarks/Alias	LDS3; mad3; LDS1C; MADH3; JV15-2; hMAD-3; hSMAD3; HSPC193; HsT17436; d3
Species	Human
Gene ID (Human)	4088
Gene ID	4088
Immunogen	Synthetic peptide corresponding to amino acids 100-200 of human Smad3 (NP_005893.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P84022
Protein Weight	48kDa
Shipping	Ice bag