


[KO Validated] Anti-Sod2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P04179

Item Number: CM0018504



Introduction

Rabbit polyclonal antibody against human SOD2, validated in WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Recognizes a 25 kDa protein. Ideal for oxidative stress research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SOD2 Rabbit pAb
Remarks/Alias	GC1, IPOB, IPO-B, MNSOD, MVCD6, GCInc1, Mn-SOD, D2
Species	Human
GeneID (human)	6648
GeneID	6648
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-222 of human SOD2 (NP_000627.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04179
Protein Weight	25 kDa
Shipping	Ice bag