

[Ko Validated] Anti-Strap Rabbit Polyclonal Antibody For Wb And Elisa - Q9Y3F4

Item Number: CM0021495

Introduction

This KO validated anti-STRAP rabbit polyclonal antibody detects STRAP in human, mouse, and rat samples via WB and ELISA. It is ideal for research on spliceosomal snRNP assembly, TGF-beta signaling, and PDPK1 regulation.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] STRAP Rabbit pAb
Remarks/Alias	MAWD; PT-WD; UNRIP; AP
Species	Human
GenelD (Human)	11171
GenelD	11171
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-350 of human STRAP (NP_009109.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y3F4
Protein Weight	38kDa
Shipping	Ice bag