

[KO Validated] Anti-Sugp2 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Ix01

Item Number: CM0022088



Introduction

Rabbit polyclonal antibody targeting human SUGP2 (SFRS14) validated for Western blot and ELISA. Detects human and mouse SUGP2, with a predicted molecular weight of 120 kDa. Suitable for mRNA splicing research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SUGP2 Rabbit pAb
Remarks/Alias	SFRS14; SRFS14; P2
Species	Human
Gene ID (Human)	10147
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 500-700 of human SUGP2 (NP_055699.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8IX01
Protein Weight	120kDa
Shipping	Ice bag