

[Ko Validated] Anti-Tfap4 Polyclonal Antibody For Wb, If/Icc, Elisa - Q01664

Item Number: CM0020775

Introduction

Rabbit polyclonal antibody targeting TFAP4 (AP-4, bHLHc41). Validated for WB, IF/ICC, ELISA. Reacts with human and rat. Ideal for studying transcription regulation via CAGCTG motif binding. For research use only.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] TFAP4 Rabbit pAb
Remarks/Alias	AP-4; bHLHc41; P4
Species	Human
GeneID (Human)	7023
GeneID	7023
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200-338 of human TFAP4 (NP_003214.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q01664
Protein Weight	39kDa
Shipping	Ice bag