

## [Ko Validated] Anti-Tmed10 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P49755

Item Number: CM0020727

### Introduction

Rabbit polyclonal antibody targeting TMED10 (TMP21/p23), validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with human, mouse, rat. Ideal for studying ER-Golgi vesicular trafficking and cargo receptor functions.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KO Validated] TMED10 Rabbit pAb                                                                                    |
| Remarks/Alias    | p23; TMP21; p24d1; S31I125; Tmp-21-I; S31III125; P24(DELTA); 10                                                     |
| Species          | Human                                                                                                               |
| Gene ID (Human)  | 10972                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 35-190 of human TMED10 (NP_006818.3). |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                            |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | P49755                                                                                                              |
| Protein Weight   | 25kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |