


[Ko Validated] Anti-Trim21/Ss-A Rabbit Monoclonal Antibody For Wb And Elisa - P19474

Item Number: CM0005836



Introduction

Rabbit monoclonal antibody against human TRIM21/SS-A (Ro52, SSA1), validated for Western blot and ELISA with cross-reactivity to rat. Suitable for studying E3 ubiquitin ligase function, autophagy, and autoimmune responses.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] TRIM21/SS-A Rabbit mAb
Remarks/Alias	SSA; RO52; SSA1; RNF81; Ro/SSA; -A
Species	Human
Gene ID (Human)	6737
Gene ID	6737
Immunogen	Recombinant protein corresponding to a sequence within amino acids 86-334 of human TRIM21 (NP_003361.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P19474
Protein Weight	54 kDa
Shipping	Ice bag