


[KO Validated] Anti-Trim25 Rabbit Monoclonal Antibody For Wb, Elisa - Q14258

Item Number: CM0012371



Introduction

Rabbit monoclonal antibody against human TRIM25 (EFP, RNF147), KO validated for Western blot and ELISA. Targets E3 ubiquitin/ISG15 ligase involved in innate immunity and antiviral defense. - Q14258

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] TRIM25 Rabbit mAb
Remarks/Alias	EFP; Z147; RNF147; ZNF147; 25
Species	Human
GeneID	7706
Immunogen	Synthetic peptide corresponding to amino acids 100-200 of human TRIM25 (Q14258)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q14258
Protein Weight	71kDa
Shipping	Ice bag