


[KO Validated] Anti-Ttc11/Fis1 Rabbit Monoclonal Antibody For Wb - Q9Y3D6

Item Number: CM0002648



Introduction

KO-validated rabbit monoclonal antibody against TTC11/FIS1, suitable for WB, IHC-P, IF/ICC, IP, and ELISA. Cross-reacts with human, mouse, and rat. Used in mitochondrial fission and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] TTC11/FIS1 Rabbit mAb
Remarks/Alias	TTC11; CGI-135; S1
Species	Human
GeneID (Human)	51024
GeneID	51024
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-122 of human TTC11/FIS1 (NP_057152.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y3D6
Protein Weight	17 kDa
Shipping	Ice bag