


[KO Validated] Anti-Ttc11/Fis1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Y3D6

Item Number: CM0019721



Introduction

Rabbit polyclonal antibody against human TTC11/FIS1, validated for WB and ELISA. Reacts with human, mouse, and rat. Suitable for mitochondrial fission, apoptosis, and peroxisomal fission studies. Gene ID: 51024.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] TTC11/FIS1 Rabbit pAb
Remarks/Alias	TTC11; CGI-135; S1
Species	Human
GeneID (Human)	51024
GeneID	51024
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human,Mouse,Rat
SWISS	Q9Y3D6
Protein Weight	17kDa
Shipping	Ice bag