


[KO Validated] Anti-Vat1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q99536

Item Number: CM0020639



Introduction

KO validated rabbit polyclonal antibody against VAT1, suitable for WB, IF/ICC, ELISA in human, mouse, and rat. Ideal for studying quinone oxidoreductase function, lipid transport, and mitochondrial dynamics.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] VAT1 Rabbit pAb
Remarks/Alias	VAT1; T1
Species	Human
Gene ID (Human)	10493
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human VAT1 (NP_006364.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99536
Protein Weight	42kDa
Shipping	Ice bag