

[KO Validated] Anti-Biii-Tubulin Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q13509

Item Number: CM0024272



Introduction

Rabbit polyclonal antibody against β III-Tubulin (TUBB3), validated for WB, IF-P, ELISA. Recognizes human, mouse, rat. Ideal for studying axon guidance, microtubule dynamics, and tumor expression.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] β III-Tubulin Rabbit pAb
Remarks/Alias	CDCBM; FEOM3; TUBB4; CDCBM1; CFEOM3; beta-4; CFEOM3A; in
Species	Human
GeneID (Human)	10381
GeneID	10381
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-260 of human β III-Tubulin/ β 3-Tubulin. (NP_006077.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13509
Protein Weight	50kDa
Shipping	Ice bag