


Ko Validated Anti-Von Hippel Lindau/Vhl Rabbit Monoclonal Antibody For Wb, Elisa - P40337

Item Number: CM0006749



Introduction

KO validated rabbit monoclonal antibody against human VHL, suitable for WB and ELISA. Detects human, mouse, rat VHL. Ideal for studying ubiquitination and tumor suppression pathways.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Von Hippel Lindau/VHL Rabbit mAb
Remarks/Alias	RCA1; VHL1; pVHL; HRCA1; Von Hippel Lindau/VHL
Species	Human
GeneID	7428
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of human Von Hippel Lindau/VHL (NP_000542.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40337
Protein Weight	24kDa
Shipping	Ice bag