


KO Validated Anti-Vrk1 Rabbit Polyclonal Antibody For Wb, IHC-P, Elisa - Q99986

Item Number: CM0025294



Introduction

KO validated rabbit polyclonal antibody against human VRK1, suitable for Western blot, IHC-P, and ELISA. Detects human, mouse, and rat VRK1, a serine/threonine kinase involved in cell cycle and DNA damage response.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] VRK1 Rabbit pAb
Remarks/Alias	PCH1, PCH1A, K1
Species	Human
GeneID (Human)	7443
GeneID	7443
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 257-396 of human VRK1 (NP_003375.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99986
Protein Weight	45kDa
Shipping	Ice bag