


[KO Validated] Anxa1 Rabbit Pab - P04083

Item Number: CM0018550



Introduction

Rabbit polyclonal antibody against Annexin A1 (ANXA1), validated by KO and suitable for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Target involved in inflammation and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] ANXA1 Rabbit pAb
Remarks/Alias	ANX1; LPC1; A1
Species	Human
GeneID (human)	301
GeneID	301
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-346 of human ANXA1 (NP_000691.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04083
Protein Weight	39kDa
Shipping	Ice bag