


[Ko Validated] Cdk5 Rabbit Monoclonal Antibody - Q00535

Item Number: CM0008793



Introduction

Rabbit monoclonal antibody against human CDK5, validated for WB, IF-P, IHC-P, and ELISA. Recognizes CDK5 from human, mouse, and rat. Suitable for studying neuronal development, differentiation, and apoptosis - Q00535.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CDK5 Rabbit mAb
Remarks/Alias	LIST7; PSSALRE
Species	Human
GeneID (Human)	1020
GeneID	1020
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 92-289 of human CDK5. (NP_004926.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00535
Protein Weight	33kDa
Shipping	Ice bag