


KO Validated Cdk7 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P50613

Item Number: CM0017844



Introduction

KO validated CDK7 rabbit polyclonal antibody suitable for WB, IF/ICC, ELISA. Recognizes human, mouse, and rat CDK7. Ideal for cell cycle and transcription research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CDK7 Rabbit pAb
Remarks/Alias	CAK; CAK1; HCAK; MO15; STK1; CDKN7; p39MO15; K7
Species	Human
Gene ID (Human)	1022
Gene ID	1022
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 211-346 of human CDK7 (NP_001790.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50613
Protein Weight	39kDa
Shipping	Ice bag