

[KO Validated] Cleaved Gasdermin E (N-Terminal) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O60443

Item Number: CM0011253



Introduction

Validated rabbit monoclonal antibody recognizing the cleaved N-terminal fragment of human Gasdermin E (GSDME/DFNA5). Suitable for Western blot, IHC-P, and ELISA. Ideal for pyroptosis and tumor immunity research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Cleaved Gasdermin E (N terminal) Rabbit mAb
Remarks/Alias	DFNA5, ICERE-1
Species	Human
Gene ID (Human)	1687
Gene ID	1687
Immunogen	Synthetic peptide corresponding to residues surrounding Asp270 of human Gasdermin E
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	O60443
Protein Weight	11 kDa/35 kDa/55 kDa
Shipping	Ice bag