


[KO Validated] Erk2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P28482

Item Number: CM0015847



Introduction

KO-validated rabbit polyclonal antibody against human ERK2 (MAPK1), suitable for Western blot, IF/ICC, and ELISA. Detects human, mouse, rat ERK2. Ideal for studying MAPK/ERK signaling pathways in cell growth, differentiation, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] ERK2 Rabbit pAb
Remarks/Alias	ERK; p38; p40; p41; ERK2; ERT1; NS13; ERK-2; MAPK2; PRKM1; PRKM2; P42MAPK; p41mapk; p42-MAPK; K2
Species	Human
Gene ID (Human)	5594
Gene ID	5594
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human ERK2 (NP_002736.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28482
Protein Weight	41kDa
Shipping	Ice bag