

[KO Validated] Fatty Acid Synthase (Fasn) Rabbit Polymab® For Wb, If/Icc, Elisa - P49327

Item Number: CM0006489



Introduction

High-quality [KO Validated] Fatty Acid Synthase (FASN) Rabbit PolymAb® for Western blot (WB), immunofluorescence (IF/ICC), and ELISA.

Validated for human, mouse, and rat. Ideal for fatty acid metabolism research and SARS-CoV-2 replication studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Fatty Acid Synthase (FASN) Rabbit PolymAb®
Remarks/Alias	FAS; OA-519; SDR27X1
Species	Human
Gene ID (Human)	2194
Gene ID	2194
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2212-2511 of human Fatty Acid Synthase (FASN) (NP_004095.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49327
Protein Weight	273kDa
Shipping	Ice bag