

[Ko Validated] Gbp2 Rabbit Polyclonal Antibody For Wb And Elisa - P32456

Item Number: CM0023524



Introduction

High-quality rabbit polyclonal antibody against human GBP2, validated by KO, for Western blot and ELISA. Ideal for innate immunity research, targeting a 67kDa interferon-inducible GTPase (P32456).

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] GBP2 Rabbit pAb
Remarks/Alias	GBP2; P2
Species	Human
GeneID (human)	2634
GeneID	2634
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P32456
Protein Weight	67kDa
Shipping	Ice bag