


Ko Validated Hif1An/Fih1 Rabbit Mab For Wb, Elisa - Q9Nwt6

Item Number: CM0010962



Introduction

Rabbit monoclonal antibody against human HIF1AN/FIH1 (Q9NWT6). Validated for Western blotting and ELISA. Ideal for studying hypoxia signaling, oxygen sensing, and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] HIF1AN/FIH1 Rabbit mAb
Remarks/Alias	FIH1; H1
Species	Human
GeneID (Human)	55662
GeneID	55662
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HIF1AN/FIH1 (Q9NWT6)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NWT6
Protein Weight	41kDa
Shipping	Ice bag