

[KO Validated] Hk1 Rabbit Pab For Wb, Ihc-P, If/Icc, Ip, Elisa - P19367

Item Number: CM0018250



Introduction

Rabbit polyclonal antibody against human hexokinase-1 (HK1), KO validated for Western blot, IHC-P, IF/ICC, IP, and ELISA. Reacts with human, mouse, and rat. Ideal for glycolysis, mitochondrial function, and innate immunity studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] HK1 Rabbit pAb
Remarks/Alias	HK; HKD; HKI; HXK1; NMSR; RP79; HMSNR; HK1-ta; HK1-tb; HK1-tc; NEDVIBA; hexokinase; K1
Species	Human
GeneID (Human)	3098
GeneID	3098
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-219 of human HK1 (NP_000179.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19367
Protein Weight	102kDa
Shipping	Ice bag