


[KO Validated] Hnrnpf Rabbit Pab For Wb, If/Icc, Elisa - P52597

Item Number: CM0026907



Introduction

KO validated rabbit polyclonal antibody targeting HNRNPF, suitable for Western blot, IF/ICC, and ELISA. Reacts with human and mouse HNRNPF, a key regulator of alternative splicing and pre-mRNA processing. Ideal for gene expression studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] HNRNPF Rabbit pAb
Remarks/Alias	HNRNPF; mcs94-1; OK/SW-cl.23; PF
Species	Human
GeneID (Human)	3185
GeneID	3185
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human HNRNPF (NP_001091678.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P52597
Protein Weight	46kDa
Shipping	Ice bag