


[Ko Validated] Hsp27/Hspb1 Rabbit Pab For Wb, If/Icc, Elisa - P04792

Item Number: CM0021303



Introduction

Rabbit polyclonal antibody against HSP27/HSPB1, validated by knockout, for Western blot, IF/ICC, and ELISA. Targets human heat shock protein beta-1, a molecular chaperone involved in stress response and cytoskeletal dynamics. - P04792

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] HSP27/HSPB1 Rabbit pAb
Remarks/Alias	CMT2F; HMN2B; HSP27; HSP28; Hsp25; SRP27; HS.76067; HEL-S-102; B1
Species	Human
Gene ID (Human)	3315
Gene ID	3315
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human HSP27/HSPB1 (NP_001531.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P04792
Protein Weight	23kDa
Shipping	Ice bag