

[Ko Validated] Ifitm3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q01628

Item Number: CM0013881

Introduction

Rabbit polyclonal antibody targeting human IFITM3 (SWISS: Q01628) for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, and rat. Suitable for antiviral innate immunity and interferon signaling research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] IFITM3 Rabbit pAb
Remarks/Alias	1-8U; IP15; DSPA2b; M3
Species	Human
Gene ID (Human)	10410
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-133 of human IFITM3 (NP_066362.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01628
Protein Weight	15kDa
Shipping	Ice bag