

[KO Validated] Integrin Alpha 2 (Itga2/Cd49B) Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P17301

Item Number: CM0031011



Introduction

KO Validated Rabbit Polyclonal Antibody against human Integrin alpha 2 (ITGA2/CD49b). Validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. For cell adhesion, collagen receptor, and virus entry research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Integrin alpha 2 (ITGA2/CD49b) Rabbit pAb
Remarks/Alias	BR; GPIa; CD49B; HPA-5; VLA-2; VLAA2; b)
Species	Human
GeneID (Human)	3673
GeneID	3673
Immunogen	Synthetic peptide: a synthetic peptide corresponding to a sequence within amino acids 30-162 of human Integrin alpha 2 (ITGA2/CD49b) (NP_002194.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17301
Protein Weight	129kDa
Shipping	Ice bag