

[KO Validated] Integrin Alpha 2 (Itga2/Cd49B) Rabbit Mab - P17301

Item Number: CM0003173



Introduction

KO-validated rabbit monoclonal antibody against human Integrin alpha-2 (ITGA2/CD49b), suitable for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Key reagent for collagen receptor and viral entry research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Integrin alpha 2 (ITGA2/CD49b) Rabbit mAb
Remarks/Alias	BR; GPIa; CD49B; HPA-5; VLA-2; VLAA2; b)
Species	Human
GeneID (human)	3673
GeneID	3673
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1082-1181 of human Integrin alpha 2 (ITGA2/CD49b) (P17301)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17301
Protein Weight	129kDa
Shipping	Ice bag