

[KO Validated] Irf3 Rabbit Mab - Q14653

Item Number: CM0008103



Introduction

Rabbit monoclonal antibody against human IRF3, validated for WB, IHC-P, IF/ICC, IP, and ELISA. Reacts with human, mouse, and rat. Ideal for studying innate immune signaling and type I interferon responses.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] IRF3 Rabbit mAb
Remarks/Alias	IIAE7
Species	Human
Gene ID (Human)	3661
Gene ID	3661
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human IRF3 (NP_001184051.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14653
Protein Weight	47kDa
Shipping	Ice bag