

[KO Validated] Jak2 Rabbit Mab For Wb, Ip, If/Icc, Ihc-P, Elisa - O60674

Item Number: CM0008477



Introduction

KO-validated recombinant rabbit monoclonal antibody against human JAK2 (131 kDa). Reacts with human, mouse, rat. Suitable for WB, IP, IF/ICC, IHC-P, ELISA. Critical for cytokine and JAK/STAT signaling research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] JAK2 Rabbit mAb
Remarks/Alias	JTK10
Species	Human
Gene ID (Human)	3717
Gene ID	3717
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 933-1132 of human JAK2 (NP_004963.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60674
Protein Weight	131 KDa
Shipping	Ice bag