

[KO Validated] Lamp2 Rabbit Polymab® - P13473

Item Number: CM0005639



Introduction

High-quality rabbit monoclonal antibody against human LAMP2, KO validated. Recognizes human, mouse, rat LAMP2. Suitable for WB, IF-P, IHC-P, ELISA. Ideal for autophagy and lysosomal function studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] LAMP2 Rabbit PolymAb®
Remarks/Alias	DND; LAMPB; CD107b; LAMP-2; LGP-96; LGP110
Species	Human
GeneID (Human)	3920
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13473
Protein Weight	44 kDa/45 kDa
Shipping	Ice bag