


[KO Validated] Ncoa2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q15596

Item Number: CM0015093



Introduction

Rabbit polyclonal antibody targeting human NCOA2, validated for WB, IHC-P, and ELISA. Cross-reacts with rat. Suitable for research on nuclear receptor coactivation, glucose metabolism, and circadian rhythm.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] NCOA2 Rabbit pAb
Remarks/Alias	SRC2, TIF2, GRIP1, KAT13C, NCoA-2, bHLHe75, A2
Species	Human
Gene ID (Human)	10499
Gene ID	10499
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 550-820 of human NCOA2 (NP_006531.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q15596
Protein Weight	159kDa
Shipping	Ice bag