


## [KO Validated] Nrf2 Rabbit Monoclonal Antibody For Wb, Ip, Elisa, Chip - Q16236

Item Number: CM0009641



### Introduction

High-quality [KO Validated] NRF2 Rabbit monoclonal antibody for WB, IP, ELISA, ChIP. Specific to human NRF2 (68kDa), a key oxidative stress transcription factor. Ideal for cytoprotection, immune response, and viral research.

[Learn More](#)

| Parameter        | Value                          |
|------------------|--------------------------------|
| Product Name     | [KO Validated] NRF2 Rabbit mAb |
| Remarks/Alias    | NRF2; HEBP1; Nrf-2; IMDDHH; F2 |
| Species          | Human                          |
| Gene ID (Human)  | 4780                           |
| Gene ID          | 4780                           |
| Immunogen        | Synthetic Peptide              |
| Source           | Rabbit                         |
| Category         | Monoclonal Antibodies          |
| Application      | WB, IP, ELISA, ChIP            |
| Cross Reactivity | Human                          |
| SWISS            | Q16236                         |
| Protein Weight   | 68kDa                          |
| Shipping         | Ice bag                        |