


[Ko Validated] Nrf2 Rabbit Monoclonal Antibody - Q16236

Item Number: CM0006626



Introduction

KO validated rabbit monoclonal antibody against human NRF2 (Q16236), suitable for Western blot, ChIP, and ELISA. Reacts with human, mouse, and rat. Ideal for oxidative stress and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] NRF2 Rabbit mAb
Remarks/Alias	NRF2; HEBP1; Nrf-2; IMDDHH
Species	Human
GeneID (Human)	4780
GeneID	4780
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 243-400 of human NRF2 (NP_006155.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16236
Protein Weight	68kDa
Shipping	Ice bag