


[KO Validated] Nsd3/Whsc1L1 Rabbit Pab For Wb, Elisa - Q9Bz95

Item Number: CM0025714



Introduction

KO validated rabbit polyclonal antibody against NSD3/WHSC1L1, a histone methyltransferase, for WB and ELISA. Cross-reacts with Human, Mouse, Rat. Ideal for epigenetic transcription studies. SWISS: Q9BZ95.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] NSD3/WHSC1L1 Rabbit pAb
Remarks/Alias	KMT3F; KMT3G; WHISTLE; WHSC1L1; pp14328; L1
Species	Human
GeneID (Human)	54904
GeneID	54904
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 400-600 of human NSD3/WHSC1L1 (NP_060248.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BZ95
Protein Weight	162kDa
Shipping	Ice bag