


[KO Validated] P53 Rabbit Pab For Wb, Ihc-P, If/Icc, Ip, Elisa - P04637

Item Number: CM0017632



Introduction

KO Validated p53 Rabbit Polyclonal Antibody (pAb) validated in WB, IHC-P, IF/ICC, IP, ELISA. Reacts with human, mouse, and rat p53. Targets tumor suppressor p53, ideal for apoptosis and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] p53 Rabbit pAb
Remarks/Alias	P53; BCC7; LFS1; BMFS5; TRP53; p53
Species	Human
GeneID (Human)	7157
GeneID	7157
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-393 of human p53 (NP_000537.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04637
Protein Weight	44kDa
Shipping	Ice bag