


[KO Validated] Ppp1Cc Rabbit Pab - P36873

Item Number: CM0025194



Introduction

Rabbit polyclonal antibody against human PPP1CC (PP1-gamma catalytic subunit), validated for WB, IF/ICC, and ELISA. Recognizes human, mouse, and rat PPP1CC (37 kDa). Essential for cell division, glycogen metabolism, and MAPK signaling.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PPP1CC Rabbit pAb
Remarks/Alias	PP1C; PP-1G; PPP1G; CC
Species	Human
Gene ID (Human)	5501
Gene ID	5501
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-323 of human PPP1CC (NP_002701.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36873
Protein Weight	37kDa
Shipping	Ice bag