


[Ko Validated] Rag2 Rabbit Pab For Wb, Elisa - P55895

Item Number: CM0025672



Introduction

KO Validated RAG2 rabbit polyclonal antibody for WB and ELISA, validated in Human and Mouse. Targets the 59kDa V(D)J recombination-activating protein 2, critical for immune receptor diversity.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] RAG2 Rabbit pAb
Remarks/Alias	RAG-2; G2
Species	Human
GeneID (Human)	5897
GeneID	5897
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 351-527 of human RAG2 (NP_000527.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P55895
Protein Weight	59kDa
Shipping	Ice bag