


[KO Validated] Rpl22 Rabbit Pab - P35268

Item Number: CM0023668



Introduction

Rabbit polyclonal antibody against human RPL22 (60S ribosomal protein L22). Validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Suitable for ribosomal protein and translation research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] RPL22 Rabbit pAb
Remarks/Alias	EAP; L22; HBP15; HBP15/L22; 22
Species	Human
GeneID (Human)	6146
GeneID	6146
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-122 of human RPL22 (NP_000974.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35268
Protein Weight	15kDa
Shipping	Ice bag