


[KO Validated] Stat5B Rabbit Pab For Wb, Ihc-P, If/Icc, Ip, Elisa - P51692

Item Number: CM0019352



Introduction

Rabbit polyclonal antibody targeting human STAT5B, validated for Western blot, immunohistochemistry, immunofluorescence, immunoprecipitation, and ELISA. Reacts with human, mouse, and rat. Ideal for signal transduction studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] STAT5B Rabbit pAb
Remarks/Alias	STAT5; GHISID2; 5B
Species	Human
Gene ID (Human)	6777
Gene ID	6777
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 586-787 of human STAT5B (NP_036580.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51692
Protein Weight	90kDa
Shipping	Ice bag