


[Ko Validated] Trim21/Ss-A Rabbit Mab For Wb, Elisa - P19474

Item Number: CM0005514



Introduction

Rabbit monoclonal antibody against human TRIM21/SS-A (52 kDa Ro protein), validated for WB and ELISA. Suitable for research on autoimmunity, ubiquitination, and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] TRIM21/SS-A Rabbit mAb
Remarks/Alias	SSA; RO52; SSA1; RNF81; Ro/SSA; -A
Species	Human
GeneID (Human)	6737
GeneID	6737
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P19474
Protein Weight	54kDa
Shipping	Ice bag