

[Ko Validated] Ttl12 Rabbit Polyclonal Antibody For Wb, Elisa - Q14166

Item Number: CM0020911

Introduction

KO Validated TTLL12 Rabbit Polyclonal Antibody for WB and ELISA, detects human, mouse, and rat TTLL12 (74kDa). Ideal for studying tubulin modification, mitosis, chromosome stability, and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] TTLL12 Rabbit pAb
Remarks/Alias	dj526l14.2; 12
Species	Human
Gene ID (Human)	23170
Gene ID	23170
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 485-644 of human TTLL12 (NP_055955.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14166
Protein Weight	74kDa
Shipping	Ice bag