

[Ko Validated] Txlna Rabbit Pab - P40222

Item Number: CM0021400

Introduction

Rabbit polyclonal antibody targeting human alpha-taxilin (TXLNA) for WB and ELISA. Recombinant protein immunogen against amino acids 1-160. Validated in knockout models. Ideal for vesicle traffic and exocytosis research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] TXLNA Rabbit pAb
Remarks/Alias	IL14; TXLN; NA
Species	Human
GeneID (Human)	200081
GeneID	200081
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human TXLNA (NP_787048.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P40222
Protein Weight	62kDa
Shipping	Ice bag