

## [KO Validated] Ube2C Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - 000762

Item Number: CM0023787



### Introduction

Rabbit polyclonal antibody against human UBE2C (UbcH10), validated for knockout studies. Suitable for WB, IF/ICC, and ELISA. UBE2C catalyzes ubiquitination and promotes APC/C-mediated degradation during mitosis.

[Learn More](#)

| Parameter        | Value                                                                                                                                   |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KO Validated] UBE2C Rabbit pAb                                                                                                         |
| Remarks/Alias    | UBCH10; dj447F3.2; 2C                                                                                                                   |
| Species          | Human                                                                                                                                   |
| GeneID (Human)   | 11065                                                                                                                                   |
| GeneID           | 11065                                                                                                                                   |
| Immunogen        | Recombinant protein   Recombinant fusion protein containing a sequence corresponding to amino acids 1-179 of human UBE2C (NP_008950.1). |
| Source           | Rabbit                                                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                                                   |
| Application      | WB, IF/ICC, ELISA                                                                                                                       |
| Cross Reactivity | Human                                                                                                                                   |
| SWISS            | O00762                                                                                                                                  |
| Protein Weight   | 20kDa                                                                                                                                   |
| Shipping         | Ice bag                                                                                                                                 |